

# Records Checklist

## GUIDED IMPLANT SURGERY

### PHOTOGRAPHS



- Full smile, full face photo, including profile photo.
- If class 2 or 3, retracted photos front side, left side, and right side IN occlusion.
- To keep lips and teeth clean, capture photos prior to taking impressions.

### MASTER CAST OR STL OF THE ARCH



- 1-4 implants: Maxillary & Mandibular master cast in PVS & bite registration.
- 5+ implants: Maxillary & Mandibular master cast in PVS & bite registration-Doctor will receive SCAN APPLIANCE (please follow the CBCT protocol)
- Endentous Arch: Use existing denture as a dual scan appliance.
  1. Denture must fit perfectly to the tissue & exhibit ideal tooth position.
  2. NO soft relines. Must use either hard reline if fit is not ideal, or reline with Blue Mousse impression material for both scans. Place 6 scan markers on the denture ([www.suremark.com](http://www.suremark.com)) or place 6 random gutta percha markers in the pink area of the labial flange. (CBCT protocol).
- Full Extraction & Placement: Must capture full palate if spare denture is to be ordered Edentulous opposing arch – *Impression capturing the teeth of the opposing dentition.—a crown and bridge impression of the denture.*

### QUESTIONS?



**TC Dental Laboratory**  
1000 NE 122nd Ave,  
Portland, OR 97230  
1-800-926-5412  
[info@tcdentallab.com](mailto:info@tcdentallab.com)  
[www.tcdentallab.com](http://www.tcdentallab.com)

### LOCATIONS



Portland, Oregon  
1000 NE 122nd Ave,  
Portland, OR 97230



Seattle, Washington  
907 N 130th St,  
Seattle, WA 98133



San Diego, California  
6440 Lusk Blvd d109,  
San Diego, CA 92121

### CBCT SCANS



- MORE than 50% metal based: **\*\*SCANN APPLIANCE WILL BE REQUIRED\*\***
- Scan of Patient OPEN biting on Cotton Rolls (*Place cotton rolls between cheeks and teeth*).